

# Temperature Baths for Viscometers

## Viscoblock VTB-250



Please consult with us regarding the Temperature baths available to match your application.

## Feature

- Directly heated sample cup, non-soiling heating section for cleanliness
- Wide temperature range, 250°C maximum working temperature
- Uniform temperature distribution, high adjustment accuracy
- Fast temperature rise
- Small, light weight

## Specifications

Working temperature range	room temperature+10~+250°C
Temperature setting accuracy	±0.6%+1 digit
External dimensions	W177 X D360 X H276mm (heating section), W200 X D250 X H180mm (power unit)
Weight	approx. 15kg
Maximum power consumption	0.65kVA
Power	AC100V 50/60Hz