

Silicone Calibration Fluids for TOKI Viscometers

TOKI has also developed silicone calibration fluids which enable more efficient Calibration of our single cylinder rotational viscometer (i.e. spindle rotor Viscometers such as our Model B).



Type	Viscosity (mPa-s) 25°C	Main Application
SC100	90	M1 ロータ (M1 rotor) BM, RB80L, TV-10M, TVB-20L
SC2000	1,800	H2 ロータ (H2 rotor) BH, RB-80H, TV-10H, TVB-20H